



Project Goals

Team 3


Project Goals




Identifying The Old System Issues



Old system Data collection process



Proposed Solution





Project Software



Software Features

User Friendly UI

The interface allows for an easy navigation

01



Live-plotting

Displays the curve using the data being collected to provide some form of indication to the results.

02



Quick Calibration

A 3-step calibration process that allows calibration through any weight

03



Simple Data Exporting

The data will be saved in a .txt file that can be opened as .csv file

04



Stress-Strain curve

Utilizing the Software ability to take inputs allows for a quick generation for the values

05




Easy to improve

The software allows for more improvements and extension for some features.


06



State machine




Calibration process



Main menu

Select the calibration option
by pressing A



Tare menu

Following the instruction
simply pull the Reel
And then press A




Calibrate

Put the weight and then
enter the value in grams
then press A

Done!

The value is now saved

Live plot



Drawing The background

The plot function start by drawing the background using pre-set Axis values.


Drawing The points

- The function utilizes transformations to compare the input value to reference axis value.
- Then it compare it to the previous value in order to figure out the position of the dot
- The final step is drawing the line between this point and the previous one.

Usable on any TFT screen


The function can be used on any screen if it uses the Adafruit graphics library

Settings menu




Settings menu

Select which Axis you would like to change



X-axis



Enter the expected max value for the Distance




Y-axis

Enter the expected max value for the weight

Live plot example



Stress-Strain menus




Enter the length in mm:

A: Done Length

#: Del

Length menu

Enter the length of the sample in mm




Enter Thickness in Micro meter:

A: Done Thickness

#: Del

Thickness menu

Enter the Thickness of the sample in micrometer



Enter the Width in mm:


A: Done Width

#: Del 30

Width menu

Enter the width of the sample in mm

Saving Data with SD card



File name is :

Test_Stress_Strain_1_1_2000_0_2_34.txt

A : Done

Why A Lazer Range Sensor instead of other possiblityts?

01

Manual

The range sensor has less of a chance to have a technological error

02

Accuracy

The range sensor has a smaller deviation from the actual measured distance

03


Instilation

Less coding required compared to other options




Why Not A Motor


Stepper Motor vs Lazer Sensor




VS



Stress-Strain Curve for Latex




Stress-Strain Curve for Nitrile Gloves




Item	Part #	Cost
Time Of Flight Sensor	VL53L1X	\$15
wire for distance sensor	STEMMA QT header	\$1
Load Cell 20kg & Amplifier	Hx711	\$9.5
Data logger	HiLetgoSD card reader	\$7
Display screen	2.2" TFT Screen	\$25
Card Reader/Writer	USB/MicroSD adaptor	\$6
Fishing reel	Zebco 606	\$21
MicroSD	Verbatim 16GB	\$5
micro-controller	Arduino mega	\$24
Inputpad	4x4 Keypad	\$13
RTC module	ds1307	\$4
TOTAL		\$152

Front



Top



Crank and Pulley System

- 3 pulley mechanical advantage
 - Supports 3x the applied load
 - Allows for loads greater than 4kg to be applied (up to 6kg)
- Failsafe to prevent damaging load cell
- Crank provides greater control and consistency when taking measurements
- Entire system costs less than \$30

Upgraded Clamps

- Much greater clamping strength than original design
 - Prevents slipping when using thin samples
- Heavy duty to grip and hold larger samples
- Much quicker and easier to secure samples than original design

Contribution

Q&A
